

**Instructions for collection and
transport of specimens for food
poisoning**

Samples for food poisoning

- 1. Nasal and Hand swab for food handler**
- 2. Rectal swab and vomitus from patients**
- 3. Remnant of cooked food**
- 4. Water sample**

(1) Collection and transport of Nasal and Hand swab

- **Use sterile cotton swab**
- **Introduce one inch into nose, gently rotate against the nasal mucosa of both sides**
- **After washing hands with soap and water (not antiseptic), rub with swab stick**
- **Place it in a transport media (Stuart or Amies)**
- **Label and send to lab as soon as possible**

Remark: Only one swab must be used for collection of 10 fingers.
Only one swab must be used for collection of both nostrils.

(2) Collection and transport of Rectal swab

- **Swab introduce at least 2 inches deep within the rectum**
- **After withdrawing of swab from rectum, check faecal stain**
- **Swab is put into transport media (Normal saline can be used)**

Collection and transport of vomitus

- **Vomitus must be collected in sterile container**

(3) Collection and transport of remnant of cooked food

- **Remnant of food must be collected in sterile (clean) plastic bag separately**
- **Proper labeling**
- **Send with cold chain**

Remark: Raw materials used (before cooking) must be sent to FDA

1. Water sample

Water sources can be divided into three basic types

- (a) Water from a tap or fixed hand pump**
- (b) Water from a reservoir (lake, tank, river)**
- (c) Water from a dug well**

Sampling from a tap or pump outlet

- **Remove any attachments from tap that may cause splashing.**
- **Wipe off the dirt from outside the tap.**
- **Turn on the tap at maximum flow and let the water flow for 1-2 minutes.**
- **Sterilize it for a minute with a flame using a gas lighter or ignited cotton wool soaked in spirit.**
- **Turn on the tap and allow the water to flow at medium flow for 1-2 minutes.**
- **Open the sterile container for collecting the sample and fill by holding the bottle under the water jet. Leave a small airspace to facilitate shaking at the time of inoculation prior to analysis.**
- **Stopper the cap and label the specimen**

Sampling from a reservoir

- **Open the bottle under sterile conditions.**
- **Fill it by holding the bottle by the lower part, submerging it to a depth of about 20 cm., with mouth facing upwards. If there is a current, the bottle should face the current.**
- **Stopper the bottle and label it.**

Sampling from a dug well

- **Attach a stone to the sampling bottle with a piece of string.**
- **Tie a 20 meter length of clean string on the bottle and to a stick.**
- **Open the bottle as described above and lower into the well.**
- **Immerse the bottle completely in water without touching the sides of the well and lower it down to the bottom of the well.**
- **Pull it out when the bottle is filled.**
- **Discard a little water to provide airspace.**
- **Stopper and label the bottle.**

The water sample should be transported to the laboratory as soon as possible, preferably within one hour. If it takes more than **three hours**, it should be transported in **ice box** and should be processed within **24 hours** of collection. While sampling chlorinated water 0.5 ml of **sodium thiosulphate solution** (18 gm/L) should be added to the sampling bottles to neutralize the residual chlorine present in water.